

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
Modified: 5/29/2025 10:41:59 AM

Electrophoresis File Run Summary

Instrument Information:





Instrument Name: DE34903152 Firmware: C.01.069
Serial#: DE34903152 Type: G2938C

Assay Information:

Assay Origin Path: C:\Program Files (x86)\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\High Sensitivity DNA.xsy
Assay Class: High Sensitivity DNA Assay
Version: 1.03
Assay Comments: Copyright © 2003-2010 Agilent Technologies

Chip Information:

Chip Lot #:
Reagent Kit Lot #:
Chip Comments:

W0586**W0587****W0588****W0589**

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
Modified: 5/29/2025 10:41:59 AM

Electrophoresis File Run Summary (Chip Summary)

Sample Name	Sample Comment	Rest. Digest	Status	Observation	Result Label	Result Color
W0586		☐	✓			
W0587		☐	✓			
W0588		☐	✓			
W0589		☐	✓			
Ladder		☐	✓			

Chip Lot #**Reagent Kit Lot #****Chip Comments :**

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
Modified: 5/29/2025 10:41:59 AM

Electrophoresis Assay Details

General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 12
Minimum Visible Range [s] : 32
Maximum Visible Range [s] : 138
Start Analysis Time Range [s] : 33
End Analysis Time Range [s] : 137.5
Ladder Concentration [pg/μl] : 1950
Uses Standard Area for Ladder Fragments
Lower Marker Concentration [pg/μl] : 125
Upper Marker Concentration [pg/μl] : 75
Used Upper Marker for Quantitation
Standard Curve Fit is Point to Point
Show Data Aligned to Lower and Upper Marker

Integrator Settings

Integration Start Time [s] : 33.05
Integration End Time [s] : 137
Slope Threshold : 0.8
Height Threshold [FU] : 5
Area Threshold : 0.1
Width Threshold [s] : 0.6
Baseline Plateau [s] : 0.5

Filter Settings

Filter Width [s] : 0.5
Polynomial Order : 4


Ladder

Ladder Peak	Size	Area
1	35	160
2	50	210
3	100	208
4	150	221
5	200	242
6	300	270
7	400	305
8	500	306
9	600	336
10	700	321
11	1000	366
12	2000	413
13	3000	411
14	7000	400
15	10380	214

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
Modified: 5/29/2025 10:41:59 AM

Electropherogram Summary



Overall Results for Ladder

Noise: 0.2


Peak table for Ladder

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	46	0.00	0.0	
3	50	150.00	4,545.5	Ladder Peak
4	100	150.00	2,272.7	Ladder Peak
5	150	150.00	1,515.2	Ladder Peak
6	200	150.00	1,136.4	Ladder Peak
7	300	150.00	757.6	Ladder Peak
8	400	150.00	568.2	Ladder Peak
9	500	150.00	454.5	Ladder Peak
10	600	150.00	378.8	Ladder Peak
11	700	150.00	324.7	Ladder Peak
12	1,000	150.00	227.3	Ladder Peak
13	2,000	150.00	113.6	Ladder Peak
14	3,000	150.00	75.8	Ladder Peak
15	7,000	150.00	32.5	Ladder Peak
16	10,380	75.00	10.9	Upper Marker

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
 Modified: 5/29/2025 10:41:59 AM

Electropherogram Summary Continued ...



Overall Results for sample 1 : W0586

Number of peaks found: 9 Corr. Area 1: 685.8
 Noise: 0.2

Peak table for sample 1 : W0586

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	476	20.88	66.5	
3	518	49.21	144.1	
4	567	107.32	287.0	
5	640	198.06	468.6	
6	873	356.22	617.9	
7	909	146.46	244.2	
8	1,247	369.68	449.2	
9	1,510	158.35	158.9	
10	1,764	675.24	580.1	
11	10,380	75.00	10.9	Upper Marker


Region table for sample 1 : W0586

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	1,000	685.8	39	728	22.7	1,040.46	2,337.9	Blue

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
 Modified: 5/29/2025 10:41:59 AM

Electropherogram Summary Continued ...



Overall Results for sample 2 : W0587

Number of peaks found: 11 Corr. Area 1: 576.4
 Noise: 0.2

Peak table for sample 2 : W0587

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	475	32.02	102.2	
3	510	39.18	116.4	
4	566	95.22	254.9	
5	637	145.66	346.6	
6	890	333.93	568.4	
7	934	93.32	151.3	
8	1,037	96.20	140.5	
9	1,232	202.21	248.7	
10	1,524	186.80	185.7	
11	1,736	266.47	232.6	
12	2,234	186.34	126.4	
13	10,380	75.00	10.9	Upper Marker


Region table for sample 2 : W0587

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	1,000	576.4	42	724	23.5	873.59	1,994.3	Blue

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
 Modified: 5/29/2025 10:41:59 AM

Electropherogram Summary Continued ...



Overall Results for sample 3 : W0588

Number of peaks found: 9 Corr. Area 1: 538.8
 Noise: 0.2

Peak table for sample 3 : W0588

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	477	21.12	67.1	
3	510	35.16	104.4	
4	563	86.87	233.6	
5	638	187.10	444.1	
6	900	473.41	797.2	
7	1,061	97.59	139.4	
8	1,223	212.88	263.7	
9	1,502	221.31	223.2	
10	1,745	660.21	573.1	
11	10,380	75.00	10.9	Upper Marker


Region table for sample 3 : W0588

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	1,000	538.8	39	731	22.1	912.84	2,024.8	Blue

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
 Modified: 5/29/2025 10:41:59 AM

Electropherogram Summary Continued ...



Overall Results for sample 4 : W0589

Number of peaks found: 6 Corr. Area 1: 668.1
 Noise: 0.2

Peak table for sample 4 : W0589

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	574	150.16	396.7	
3	634	196.16	468.4	
4	897	527.28	890.6	
5	1,225	327.46	404.9	
6	1,482	251.57	257.2	
7	1,684	734.47	661.0	
8	10,380	75.00	10.9	Upper Marker

Region table for sample 4 : W0589

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	1,000	668.1	41	721	23.2	1,141.25	2,602.4	Blue

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
Modified: 5/29/2025 10:41:59 AM

Gel Image

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2025-05-29\10XSC0427 cDNA.xad

Created: 5/29/2025 9:52:00 AM
 Modified: 5/29/2025 10:41:59 AM

Run Logbook

Description	Number	Source	Category	Sub Category	Time	Time Zone	User	Host
Run ended on port 1 (Number of wells acquired: 12)		Instrument	Run		5/29/2025 10:32:27 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Run started on port 1 (File: C:\Program Files (x86)\Agilent\2100 bioanalyzer\2100 expert\data\2025-05-29\Bioanalyzer1_High Sensitivity DNA Assay_DE34903152_2025-05-29_002.xad)		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Product Number : G2938C		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Name :		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Vendor : Agilent Technologies		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Serial# : DE34903152		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Firmware : C.01.069		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Cartridge : Electrode		Instrument	Run		5/29/2025 9:52:06 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB