


Assay Class: High Sensitivity DNA Assay
Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
Modified: 6/2/2023 1:02:46 PM

Electrophoresis File Run Summary



Instrument Information:

Instrument Name: DE34903152 Firmware: C.01.069
Serial#: DE34903152 Type: G2938C

Assay Information:

Assay Origin Path: C:\Program Files (x86)\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\High Sensitivity DNA.xsy
Assay Class: High Sensitivity DNA Assay
Version: 1.03
Assay Comments: Copyright © 2003-2010 Agilent Technologies


Chip Information:

Chip Lot #:
Reagent Kit Lot #:
Chip Comments:


PB921_CMPRA3_Amplicon_QC_pre-lib prep

PB923_NIJL009_Amplicon_QC_pre-lib prep


PB929_CH1_Amplicon_QC_pre-lib prep




QC1500_CT1




QC1501_S1-2




RN010_QC




QC1497_01_LWB01




QC1497_02_LWB02




QC1497_03_LWB03




QC1497_04_LWB04



QC1497_05_LWB05



QC1497_04_LWB04



QC1497_05_LWB05

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
 Modified: 6/2/2023 1:02:46 PM

Electrophoresis File Run Summary (Chip Summary)

Sample Name	Sample Comment	Rest. Digest	Status	Observation	Result Label	Result Color
PB921_CMPRA3_Amplicon_QC_pre-lib prep		<input type="checkbox"/>	✓			
PB923_NIJL009_Amplicon_QC_pre-lib prep		<input type="checkbox"/>	✓			
PB929_CH1_Amplicon_QC_pre-lib prep		<input type="checkbox"/>	✓			
QC1500_CT1		<input type="checkbox"/>	✓			
QC1501_S1-2		<input type="checkbox"/>	✓			
RN010_QC		<input type="checkbox"/>	✓			
QC1497_01_LWB01		<input type="checkbox"/>	✓			
QC1497_02_LWB02		<input type="checkbox"/>	✓			
QC1497_03_LWB03		<input type="checkbox"/>	✓			
QC1497_04_LWB04		<input type="checkbox"/>	✓			
QC1497_05_LWB05		<input type="checkbox"/>	✓			
Ladder		<input type="checkbox"/>	✓			

Chip Lot #

Reagent Kit Lot #

Chip Comments :

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
 Modified: 6/2/2023 1:02:46 PM

Electrophoresis Assay Details

General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 12
 Minimum Visible Range [s] : 32
 Maximum Visible Range [s] : 138
 Start Analysis Time Range [s] : 33
 End Analysis Time Range [s] : 137.5
 Ladder Concentration [pg/μl] : 1950
 Uses Standard Area for Ladder Fragments
 Lower Marker Concentration [pg/μl] : 125
 Upper Marker Concentration [pg/μl] : 75
 Used Upper Marker for Quantitation
 Standard Curve Fit is Point to Point
 Show Data Aligned to Lower and Upper Marker

Integrator Settings

Integration Start Time [s] : 33.05
 Integration End Time [s] : 137
 Slope Threshold : 0.8
 Height Threshold [FU] : 5
 Area Threshold : 0.1
 Width Threshold [s] : 0.6
 Baseline Plateau [s] : 0.5

Filter Settings

Filter Width [s] : 0.5
 Polynomial Order : 4


Ladder

Ladder Peak	Size	Area
1	35	160
2	50	210
3	100	208
4	150	221
5	200	242
6	300	270
7	400	305
8	500	306
9	600	336
10	700	321
11	1000	366
12	2000	413
13	3000	411
14	7000	400
15	10380	214

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
 Modified: 6/2/2023 1:02:46 PM

Electropherogram Summary



Overall Results for Ladder

Noise: 0.2


Peak table for Ladder

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	50	150.00	4,545.5	Ladder Peak
3	100	150.00	2,272.7	Ladder Peak
4	150	150.00	1,515.2	Ladder Peak
5	200	150.00	1,136.4	Ladder Peak
6	300	150.00	757.6	Ladder Peak
7	400	150.00	568.2	Ladder Peak
8	500	150.00	454.5	Ladder Peak
9	600	150.00	378.8	Ladder Peak
10	700	150.00	324.7	Ladder Peak
11	1,000	150.00	227.3	Ladder Peak
12	2,000	150.00	113.6	Ladder Peak
13	3,000	150.00	75.8	Ladder Peak
14	7,000	150.00	32.5	Ladder Peak
15	10,380	75.00	10.9	Upper Marker

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
 Modified: 6/2/2023 1:02:46 PM

Electropherogram Summary Continued ...



Overall Results for sample 2 : PB923 NIJL009 Amplicon QC pre-lib prep

Number of peaks found: 4 Corr. Area 1: 1,147.1
 Noise: 0.2

Peak table for sample 2 : PB923 NIJL009 Amplicon QC pre-lib prep

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	2,560	4.14	2.5	excluded peak
3	3,402	30.99	13.8	
4	5,173	270.65	79.3	
5	6,202	228.27	55.8	
6	10,380	75.00	10.9	Upper Marker
7	13,895	0.00	0.0	


Region table for sample 2 : PB923 NIJL009 Amplicon QC pre-lib prep

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
2,688	9,222	1,147.1	97	5,642	16.7	529.92	143.9	Blue

Assay Class: High Sensitivity DNA Assay
Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
Modified: 6/2/2023 1:02:46 PM

Gel Image




Assay Class: High Sensitivity DNA Assay
Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
Modified: 6/2/2023 1:02:46 PM

Curves

Standard Curve



Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad

Created: 6/2/2023 11:48:05 AM
 Modified: 6/2/2023 1:02:46 PM

Run Logbook

Description	Number	Source	Category	Sub Category	Time	Time Zone	User	Host
Run ended on port 1 (Number of wells acquired: 12)		Instrument	Run		6/2/2023 12:29:28 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Run started on port 1 (File: C:\Users\sbsuser\Desktop\2023-06-02\Bioanalyzer1_High Sensitivity DNA Assay_2023-06-02_001.xad)		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Product Number : G2938C		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Name :		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Vendor : Agilent Technologies		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Serial# : DE34903152		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Firmware : C.01.069		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Cartridge : Electrode		Instrument	Run		6/2/2023 11:48:11 AM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB