


Assay Class: High Sensitivity DNA Assay  
Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
Modified: 5/4/2023 2:52:10 PM

**Electrophoresis File Run Summary**



Instrument Information:






Instrument Name: DE34903152      Firmware: C.01.069  
Serial#: DE34903152      Type: G2938C

Assay Information:

Assay Origin Path: C:\Program Files (x86)\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\High Sensitivity DNA.xsy  
Assay Class: High Sensitivity DNA Assay  
Version: 1.03  
Assay Comments: Copyright © 2003-2010 Agilent Technologies

Chip Information:

Chip Lot #:  
Reagent Kit Lot #:  
Chip Comments:



Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electrophoresis File Run Summary (Chip Summary)**

Sample Name	Sample Comment	Rest. Digest	Status	Observation	Result Label	Result Color
MM48	QC1467	<input type="checkbox"/>	✓			
MM72	QC1467	<input type="checkbox"/>	✓			
MB40	QC1467	<input type="checkbox"/>	✓			
MB80	QC1467	<input type="checkbox"/>	✓			
MB108	QC1467	<input type="checkbox"/>	✓			
Ladder		<input type="checkbox"/>	✓			

**Chip Lot #**

**Reagent Kit Lot #**

**Chip Comments :**

Assay Class: High Sensitivity DNA Assay  
Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
Modified: 5/4/2023 2:52:10 PM

## Electrophoresis Assay Details

### General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 12  
Minimum Visible Range [s] : 32  
Maximum Visible Range [s] : 138  
Start Analysis Time Range [s] : 33  
End Analysis Time Range [s] : 137.5  
Ladder Concentration [pg/ $\mu$ l] : 1950  
Uses Standard Area for Ladder Fragments  
Lower Marker Concentration [pg/ $\mu$ l] : 125  
Upper Marker Concentration [pg/ $\mu$ l] : 75  
Used Upper Marker for Quantitation  
Standard Curve Fit is Point to Point  
Show Data Aligned to Lower and Upper Marker

### Integrator Settings

Integration Start Time [s] : 33.05  
Integration End Time [s] : 137  
Slope Threshold : 0.8  
Height Threshold [FU] : 5  
Area Threshold : 0.1  
Width Threshold [s] : 0.6  
Baseline Plateau [s] : 0.5

### Filter Settings

Filter Width [s] : 0.5  
Polynomial Order : 4


### Ladder

Ladder Peak	Size	Area
1	35	160
2	50	210
3	100	208
4	150	221
5	200	242
6	300	270
7	400	305
8	500	306
9	600	336
10	700	321
11	1000	366
12	2000	413
13	3000	411
14	7000	400
15	10380	214

Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electropherogram Summary**



**Overall Results for Ladder**

Noise: 0.3


**Peak table for Ladder**

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	50	150.00	4,545.5	Ladder Peak
3	100	150.00	2,272.7	Ladder Peak
4	150	150.00	1,515.2	Ladder Peak
5	200	150.00	1,136.4	Ladder Peak
6	300	150.00	757.6	Ladder Peak
7	400	150.00	568.2	Ladder Peak
8	500	150.00	454.5	Ladder Peak
9	600	150.00	378.8	Ladder Peak
10	700	150.00	324.7	Ladder Peak
11	1,000	150.00	227.3	Ladder Peak
12	2,000	150.00	113.6	Ladder Peak
13	3,000	150.00	75.8	Ladder Peak
14	7,000	150.00	32.5	Ladder Peak
15	10,380	75.00	10.9	Upper Marker

Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electropherogram Summary Continued ...**



**Overall Results for sample 7 : MM48**

Number of peaks found: 1                      Corr. Area 1: 1,044.9  
 Noise: 0.3

**Peak table for sample 7 : MM48**

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	619	1,062.83	2,600.2	
3	10,380	75.00	10.9	Upper Marker


**Region table for sample 7 : MM48**

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
410	1,740	1,044.9	85	643	16.9	1,089.42	2,614.4	Blue

Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electropherogram Summary Continued ...**



**Overall Results for sample 8 : MM72**

Number of peaks found: 18                      Corr. Area 1: 1,037.9  
 Noise: 0.4

**Peak table for sample 8 : MM72**

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	179	9.44	79.8	
3	195	63.47	493.6	
4	204	63.13	469.8	
5	235	52.81	340.8	
6	259	31.89	186.6	
7	275	29.88	164.4	
8	292	21.41	111.1	
9	312	8.11	39.4	
10	327	7.50	34.8	
11	366	9.23	38.2	
12	384	8.89	35.1	
13	408	14.17	52.6	
14	434	13.35	46.6	
15	461	32.29	106.2	
16	493	9.96	30.6	
17	623	881.27	2,144.6	
18	2,535	2.81	1.7	
19	10,380	75.00	10.9	Upper Marker
20	13,164	0.00	0.0	


**Region table for sample 8 : MM72**

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
399	1,963	1,037.9	69	661	29.5	1,007.88	2,449.3	Blue

Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electropherogram Summary Continued ...**



**Overall Results for sample 9 : MB40**

Number of peaks found: 4                      Corr. Area 1: 851.7  
 Noise: 0.2

**Peak table for sample 9 : MB40**

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	424	254.68	909.2	
3	446	692.91	2,353.0	
4	543	25.53	71.3	
5	601	9.33	23.5	
6	10,380	75.00	10.9	Upper Marker


**Region table for sample 9 : MB40**

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
293	1,000	851.7	74	464	17.3	1,105.77	3,699.5	Blue

Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electropherogram Summary Continued ...**



**Overall Results for sample 10 : MB80**

Number of peaks found: 10      Corr. Area 1: 829.1  
 Noise: 0.3

**Peak table for sample 10 : MB80**

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	444	645.58	2,204.0	
3	540	9.03	25.3	
4	555	12.70	34.7	
5	606	11.04	27.6	
6	2,348	17.94	11.6	
7	2,624	5.68	3.3	
8	2,901	8.36	4.4	
9	3,628	8.10	3.4	
10	8,027	8.58	1.6	
11	10,380	75.00	10.9	Upper Marker
12	12,793	0.00	0.0	

**Region table for sample 10 : MB80**


From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
299	760	829.1	77	455	11.3	704.99	2,375.7	Blue



Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Electropherogram Summary Continued ...**



**Overall Results for sample 11 : MB108**

Number of peaks found: 8      Corr. Area 1: 850.6  
 Noise: 0.3

**Peak table for sample 11 : MB108**

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	441	645.93	2,216.9	
3	540	16.51	46.4	
4	577	7.67	20.1	
5	603	18.54	46.6	
6	2,342	29.42	19.0	
7	4,141	15.62	5.7	
8	6,329	9.40	2.2	
9	8,130	8.17	1.5	
10	10,380	75.00	10.9	Upper Marker


**Region table for sample 11 : MB108**

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
299	1,475	850.6	79	481	29.4	743.73	2,464.3	Blue

Assay Class: High Sensitivity DNA Assay  
Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
Modified: 5/4/2023 2:52:10 PM

**Gel Image**



Assay Class: High Sensitivity DNA Assay  
 Data Path: C:\Users\sbsuser\Desktop\2023-05-04\QC1467\_BioA.xad

Created: 5/4/2023 1:15:52 PM  
 Modified: 5/4/2023 2:52:10 PM

**Run Logbook**

Description	Number	Source	Category	Sub Category	Time	Time Zone	User	Host
Run ended on port 1 (Number of wells acquired: 12)		Instrument	Run		5/4/2023 1:56:14 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Run started on port 1 (File: C:\Users\sbsuser\Desktop\2023-05-04\Bioanalyzer1_High Sensitivity DNA Assay_2023-05-04_002.xad)		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Product Number : G2938C		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Name :		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Vendor : Agilent Technologies		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Serial# : DE34903152		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Firmware : C.01.069		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB
Cartridge : Electrode		Instrument	Run		5/4/2023 1:15:53 PM	(GMT --07:00) Pacific Standard Time	sbsuser	DESKTOP-4UNV VOB