


Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electrophoresis File Run Summary



Instrument Information:

Instrument Name: DE13701086 Firmware: C.01.069
Serial#: DE13701086 Type: G2938B

Assay Information:


Assay Origin Path: C:\Program Files\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\DNA 12000 Series II.xsy

Assay Class: DNA 12000
Version: 2.4
Assay Comments: DNA Analysis 100 -12000 bp

© Copyright 2003-2009 Agilent Technologies, Inc.

Chip Information:

Chip Lot #: Reagent Kit Lot #: Chip Comments:



Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electrophoresis File Run Summary (Chip Summary)

Sample Name	Sample Comment	Rest. Digest	Status	Observation	Result Label	Result Color
NC #5	PacBio Workshop	<input type="checkbox"/>	✓			
Sample 2		<input type="checkbox"/>	✓			
JZ	PacBio Workshop	<input type="checkbox"/>	✓			
Sample 4		<input type="checkbox"/>	✓			
IH	PacBio Workshop	<input type="checkbox"/>	✓			
Sample 6		<input type="checkbox"/>	✓			
VL	PacBio Workshop	<input type="checkbox"/>	✓			
Sample 8		<input type="checkbox"/>	✓			
SCL	PacBio Workshop	<input type="checkbox"/>	✓			
Sample 10		<input type="checkbox"/>	✓			
#2	PacBio Workshop	<input type="checkbox"/>	✓			
Sample 12		<input type="checkbox"/>	✓			
Ladder		<input type="checkbox"/>	✓			

Chip Lot #

Reagent Kit Lot #

Chip Comments :

Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electrophoresis Assay Details

General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 13
 Minimum Visible Range [s] : 20
 Maximum Visible Range [s] : 99
 Start Analysis Time Range [s] : 20
 End Analysis Time Range [s] : 98.95
 Ladder Concentration [ng/μl] : 44
 Uses Standard Area for Ladder Fragments
 Lower Marker Concentration [ng/μl] : 8.3
 Upper Marker Concentration [ng/μl] : 4.2
 Used Upper Marker for Quantitation
 This is a Qualitative Assay Only
 Standard Curve Fit is Point to Point
 Show Data Aligned to Lower and Upper Marker

Integrator Settings

Integration Start Time [s] : 20
 Integration End Time [s] : 98.95
 Slope Threshold : 0.8
 Height Threshold [FU] : 20
 Area Threshold : 0.1
 Width Threshold [s] : 0.5
 Baseline Plateau [s] : 0.5

Filter Settings

Filter Width [s] : 0.5
 Polynomial Order : 4


Ladder

Ladder Peak	Size	Area
1	50	120
2	100	48
3	300	61
4	500	77
5	700	81
6	1000	86
7	1500	92
8	2000	92
9	3000	96
10	5000	100
11	7000	98
12	10380	102
13	17000	110

Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary



Peak table for Ladder


Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	100	4.00	60.6	Ladder Peak
3	300	4.00	20.2	Ladder Peak
4	500	4.00	12.1	Ladder Peak
5	700	4.00	8.7	Ladder Peak
6	1,000	4.00	6.1	Ladder Peak
7	1,500	4.00	4.0	Ladder Peak
8	2,000	4.00	3.0	Ladder Peak
9	3,000	4.00	2.0	Ladder Peak
10	5,000	4.00	1.2	Ladder Peak
11	7,000	4.00	0.9	Ladder Peak
12	10,380	4.00	0.6	Ladder Peak
13	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

NC #5 [PacBio Workshop]



Overall Results for sample 1 : NC #5

Number of peaks found: 1

Peak table for sample 1 : NC #5


Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	11,046	5.79	0.8	
3	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

Sample 2



Overall Results for sample 2 : Sample 2

Number of peaks found: 0

Peak table for sample 2 : Sample 2


Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

JZ [PacBio Workshop]



Overall Results for sample 3 : JZ

Number of peaks found: 1

Peak table for sample 3 : JZ


Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	6,271	20.56	5.0	
3	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

Sample 4



Overall Results for sample 4 : Sample 4

Number of peaks found: 0


Peak table for sample 4 : Sample 4

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...



Overall Results for sample 5 : IH

Number of peaks found: 0 Area 1: 0.0

Peak table for sample 5 : IH

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker
3	70,645	0.00	0.0	

Region table for sample 5 : IH


From [bp]	To [bp]	Average Size [bp]	Conc. [ng/μl]	Area	% of Total	Size distribution in CV [%]	Color
2,949	32,119	13,922	0.00	0.0	0	46.5	Blue

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

Sample 6



Overall Results for sample 6 : Sample 6

Number of peaks found: 0


Peak table for sample 6 : Sample 6

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...



Overall Results for sample 7 : VL

Number of peaks found: 2

Peak table for sample 7 : VL


Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	6,880	16.55	3.6	
3	8,645	9.16	1.6	
4	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

Sample 8



Setpoint Deviations for sample 8 : Sample 8

Height Threshold [FU] : 5

Overall Results for sample 8 : Sample 8

Number of peaks found: 0


Peak table for sample 8 : Sample 8

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...



Overall Results for sample 9 : SCL

Number of peaks found: 0 Area 1: 91.8

Peak table for sample 9 : SCL

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Region table for sample 9 : SCL



From [bp]	To [bp]	Average Size [bp]	Conc. [ng/μl]	Area	% of Total	Size distribution in CV [%]	Color
1,338	36,130	13,053	7.98	91.8	78	59.1	Blue

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

Sample 10



Setpoint Deviations for sample 10 : Sample 10

Height Threshold [FU] : 5

Overall Results for sample 10 : Sample 10

Number of peaks found: 0


Peak table for sample 10 : Sample 10

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...



Overall Results for sample 11 : #2

Number of peaks found: 0 Area 1: 78.5

Peak table for sample 11 : #2

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Region table for sample 11 : #2


From [bp]	To [bp]	Average Size [bp]	Conc. [ng/μl]	Area	% of Total	Size distribution in CV [%]	Color
1,164	46,075	12,987	1.88	78.5	68	65.3	Blue

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Electropherogram Summary Continued ...

Sample 12



Setpoint Deviations for sample 12 : Sample 12

Height Threshold [FU] : 5

Overall Results for sample 12 : Sample 12

Number of peaks found: 0


Peak table for sample 12 : Sample 12

Peak	Size [bp]	Conc. [ng/μl]	Molarity [nmol/l]	Observations
1	50	8.30	251.5	Lower Marker
2	17,000	4.20	0.4	Upper Marker

Assay Class: DNA 12000
Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
Modified: 4/14/2015 2:44:28 PM

Gel Image



Assay Class: DNA 12000
 Data Path: C:\...ents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad

Created: 4/14/2015 1:49:00 PM
 Modified: 4/14/2015 2:44:28 PM

Run Logbook

Description	Number	Source	Category	Sub Category	Time	Time Zone	User	Host
Run ended on port 1 (Number of wells acquired: 13)		Instrument	Run		4/14/2015 2:24:01 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Run started on port 1 (File: C:\Documents and Settings\Bioanalyzer\2015-04-14\2015-04-14_002.xad)		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Product Number : G2938B		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Name :		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Vendor : Agilent Technologies		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Serial# : DE13701086		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Firmware : C.01.069		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1
Cartridge : Electrode		Instrument	Run		4/14/2015 1:49:05 PM	(GMT --07:00) Pacific Standard Time	UC Davis	D8XSMGH1