

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
Modified: 7/10/2017 1:02:53 PM

Electrophoresis File Run Summary

Instrument Information:


Instrument Name: DE13701086 Firmware: C.01.069
Serial#: DE13701086 Type: G2938B

Assay Information:

Assay Origin Path: C:\Program Files\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\High Sensitivity DNA.xsy
Assay Class: High Sensitivity DNA Assay
Version: 1.03
Assay Comments: Copyright © 2003-2010 Agilent Technologies

Chip Information:

Chip Lot #:
Reagent Kit Lot #:
Chip Comments:



Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electrophoresis File Run Summary (Chip Summary)

Sample Name	Sample Comment	Rest. Digest	Status	Observation	Result Label	Result Color
2039-2		<input type="checkbox"/>	✓			
2040-2		<input type="checkbox"/>	✓			
2048-2		<input type="checkbox"/>	✓			
2050-2		<input type="checkbox"/>	✓			
2053-2		<input type="checkbox"/>	✓			
2055-2		<input type="checkbox"/>	✓			
2056-2		<input type="checkbox"/>	✓			
2057-2		<input type="checkbox"/>	✓			
2058-2		<input type="checkbox"/>	✓			
2063-2		<input type="checkbox"/>	✓			
2069-2		<input type="checkbox"/>	✓			
Ladder		<input type="checkbox"/>	✓			

Chip Lot #

Reagent Kit Lot #

Chip Comments :

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
Modified: 7/10/2017 1:02:53 PM

Electrophoresis Assay Details

General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 12
Minimum Visible Range [s] : 32
Maximum Visible Range [s] : 138
Start Analysis Time Range [s] : 33
End Analysis Time Range [s] : 137.5
Ladder Concentration [pg/μl] : 1950
Uses Standard Area for Ladder Fragments
Lower Marker Concentration [pg/μl] : 125
Upper Marker Concentration [pg/μl] : 75
Used Upper Marker for Quantitation
Standard Curve Fit is Point to Point
Show Data Aligned to Lower and Upper Marker

Integrator Settings

Integration Start Time [s] : 33.05
Integration End Time [s] : 137
Slope Threshold : 0.8
Height Threshold [FU] : 5
Area Threshold : 0.1
Width Threshold [s] : 0.6
Baseline Plateau [s] : 0.5

Filter Settings

Filter Width [s] : 0.5
Polynomial Order : 4


Ladder

Ladder Peak	Size	Area
1	35	160
2	50	210
3	100	208
4	150	221
5	200	242
6	300	270
7	400	305
8	500	306
9	600	336
10	700	321
11	1000	366
12	2000	413
13	3000	411
14	7000	400
15	10380	214

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary



Overall Results for Ladder

Noise: 0.1


Peak table for Ladder

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	50	150.00	4,545.5	Ladder Peak	45.33
3	100	150.00	2,272.7	Ladder Peak	50.92
4	150	150.00	1,515.2	Ladder Peak	55.66
5	200	150.00	1,136.4	Ladder Peak	60.36
6	300	150.00	757.6	Ladder Peak	69.52
7	400	150.00	568.2	Ladder Peak	77.58
8	500	150.00	454.5	Ladder Peak	83.27
9	600	150.00	378.8	Ladder Peak	87.82
10	700	150.00	324.7	Ladder Peak	90.94
11	1,000	150.00	227.3	Ladder Peak	94.94
12	2,000	150.00	113.6	Ladder Peak	101.37
13	3,000	150.00	75.8	Ladder Peak	104.59
14	7,000	150.00	32.5	Ladder Peak	109.78
15	10,380	75.00	10.9	Upper Marker	113.00

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 1 : 2039-2

Number of peaks found: 2 Corr. Area 1: 31.5
 Noise: 0.2

Peak table for sample 1 : 2039-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	153	2,163.45	21,487.9		55.90
3	174	2,240.25	19,535.4		57.90
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 1 : 2039-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	242	31.5	667.4	99.41	2	37.8

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 2 : 2040-2

Number of peaks found: 2 Corr. Area 1: 13.6
 Noise: 0.1

Peak table for sample 2 : 2040-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	151	611.72	6,121.9		55.80
3	172	778.44	6,869.9		57.70
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 2 : 2040-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	267	13.6	192.7	28.62	2	58.3

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 3 : 2048-2

Number of peaks found: 2 Corr. Area 1: 52.4
 Noise: 0.1

Peak table for sample 3 : 2048-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	152	572.53	5,699.8		55.87
3	173	1,677.05	14,708.0		57.80
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 3 : 2048-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	304	52.4	823.4	138.56	7	45.4

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 4 : 2050-2

Number of peaks found: 2 Corr. Area 1: 39.8
 Noise: 0.1

Peak table for sample 4 : 2050-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	152	1,567.82	15,601.1		55.88
3	172	2,168.99	19,107.5		57.73
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 4 : 2050-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	273	39.8	703.4	112.53	3	40.0

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 5 : 2053-2

Number of peaks found: 2 Corr. Area 1: 52.8
 Noise: 0.1

Peak table for sample 5 : 2053-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	153	861.04	8,544.4		55.92
3	174	1,566.29	13,659.2		57.90
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 5 : 2053-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	324	52.8	869.6	141.05	6	58.0

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 6 : 2055-2

Number of peaks found: 2 Corr. Area 1: 33.3
 Noise: 0.1

Peak table for sample 6 : 2055-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	151	1,609.07	16,110.2		55.79
3	172	1,307.14	11,511.3		57.74
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 6 : 2055-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	342	33.3	516.3	89.05	3	51.7

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 7 : 2056-2

Number of peaks found: 2 Corr. Area 1: 23.5
 Noise: 0.1

Peak table for sample 7 : 2056-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	152	552.64	5,517.8		55.83
3	172	1,472.28	13,001.6		57.69
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 7 : 2056-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	312	23.5	451.7	72.03	4	55.3

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 8 : 2057-2

Number of peaks found: 2 Corr. Area 1: 14.0
 Noise: 0.1

Peak table for sample 8 : 2057-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	152	1,236.03	12,335.6		55.84
3	172	1,307.17	11,521.6		57.72
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 8 : 2057-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	212	14.0	340.5	47.23	2	15.3

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 9 : 2058-2

Number of peaks found: 2 Corr. Area 1: 8.7
 Noise: 0.1

Peak table for sample 9 : 2058-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	148	705.91	7,205.7		55.52
3	170	604.36	5,376.3		57.57
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 9 : 2058-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	207	8.7	157.0	21.50	2	3.9

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 10 : 2063-2

Number of peaks found: 2 Corr. Area 1: 17.9
 Noise: 0.1

Peak table for sample 10 : 2063-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	150	1,580.89	16,003.9		55.63
3	170	1,456.25	12,976.8		57.55
4	10,380	75.00	10.9	Upper Marker	113.00


Region table for sample 10 : 2063-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	306	17.9	405.6	62.63	2	61.9

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
 Modified: 7/10/2017 1:02:53 PM

Electropherogram Summary Continued ...



Overall Results for sample 11 : 2069-2

Number of peaks found: 2 Corr. Area 1: 14.3
 Noise: 0.1

Peak table for sample 11 : 2069-2

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations	Aligned Migration Time [s]
1	35	125.00	5,411.3	Lower Marker	43.00
2	147	1,303.61	13,405.8		55.41
3	168	1,727.21	15,537.0		57.40
4	10,380	75.00	10.9	Upper Marker	113.00

Region table for sample 11 : 2069-2

From [bp]	To [bp]	Average Size [bp]	Corr. Area	Molarity [pmol/l]	Co Conc. lor [pg/μl]	% of Total	Size distribution in CV [%]
200	1,000	220	14.3	469.6	67.13	2	15.0

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
Modified: 7/10/2017 1:02:53 PM

Gel Image

Assay Class: High Sensitivity DNA Assay
Data Path: C:\... bioanalyzer\2100 expert\data\2017-07-10\2017-07-10_001.xad

Created: 7/10/2017 12:21:34 PM
Modified: 7/10/2017 1:02:53 PM

Curves

Standard Curve

