


Assay Class: High Sensitivity DNA Assay
Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
Modified: 2/18/2014 12:01:55 PM

Electrophoresis File Run Summary



Instrument Information:


Instrument Name: DE13701086 Firmware: C.01.069
Serial#: DE13701086 Type: G2938B

Assay Information:

Assay Origin Path: C:\Program Files\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\High Sensitivity DNA.xsy
Assay Class: High Sensitivity DNA Assay
Version: 1.03
Assay Comments: Copyright © 2003-2010 Agilent Technologies

Chip Information:

Chip Lot #: Reagent Kit Lot #: Chip Comments:



Assay Class: High Sensitivity DNA Assay
Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
Modified: 2/18/2014 12:01:55 PM

Electrophoresis File Run Summary (Chip Summary)

Sample Name	Sample Comment	Rest. Digest	Status	Observation	Result Label	Result Color
sample 3		<input type="checkbox"/>	✓			
sample 6		<input type="checkbox"/>	✓			
sample 11		<input type="checkbox"/>	✓			
sample 13		<input type="checkbox"/>	✓			
sample 14		<input type="checkbox"/>	✓			
sample 15		<input type="checkbox"/>	✓			
sample 17		<input type="checkbox"/>	✓			
sample 18		<input type="checkbox"/>	✓			
sample 9		<input type="checkbox"/>				
sample 10		<input type="checkbox"/>				
sample 11		<input type="checkbox"/>				
Ladder		<input type="checkbox"/>	✓			

Chip Lot #

Reagent Kit Lot #

Chip Comments :

Assay Class: High Sensitivity DNA Assay
Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
Modified: 2/18/2014 12:01:55 PM

Electrophoresis Assay Details

General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 12
Minimum Visible Range [s] : 32
Maximum Visible Range [s] : 138
Start Analysis Time Range [s] : 33
End Analysis Time Range [s] : 137.5
Ladder Concentration [pg/μl] : 1950
Uses Standard Area for Ladder Fragments
Lower Marker Concentration [pg/μl] : 125
Upper Marker Concentration [pg/μl] : 75
Used Upper Marker for Quantitation
Standard Curve Fit is Point to Point
Show Data Aligned to Lower and Upper Marker

Integrator Settings

Integration Start Time [s] : 33.05
Integration End Time [s] : 137
Slope Threshold : 0.8
Height Threshold [FU] : 5
Area Threshold : 0.1
Width Threshold [s] : 0.6
Baseline Plateau [s] : 0.5

Filter Settings

Filter Width [s] : 0.5
Polynomial Order : 4


Ladder

Ladder Peak	Size	Area
1	35	160
2	50	210
3	100	208
4	150	221
5	200	242
6	300	270
7	400	305
8	500	306
9	600	336
10	700	321
11	1000	366
12	2000	413
13	3000	411
14	7000	400
15	10380	214

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary



Overall Results for Ladder

Noise: 0.2


Peak table for Ladder

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	45.33	50	150.00	4,545.5	Ladder Peak
3	50.75	100	150.00	2,272.7	Ladder Peak
4	55.37	150	150.00	1,515.2	Ladder Peak
5	59.98	200	150.00	1,136.4	Ladder Peak
6	69.04	300	150.00	757.6	Ladder Peak
7	77.09	400	150.00	568.2	Ladder Peak
8	82.76	500	150.00	454.5	Ladder Peak
9	87.51	600	150.00	378.8	Ladder Peak
10	90.73	700	150.00	324.7	Ladder Peak
11	94.71	1,000	150.00	227.3	Ladder Peak
12	101.27	2,000	150.00	113.6	Ladder Peak
13	104.49	3,000	150.00	75.8	Ladder Peak
14	107.79	5,600	0.00	0.0	
15	109.57	7,000	150.00	32.5	Ladder Peak
16	113.00	10,380	75.00	10.9	Upper Marker
17	114.69	12,049	0.00	0.0	
18	115.46	12,800	0.00	0.0	
19	116.05	13,384	0.00	0.0	

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 1 : sample 3

Number of peaks found: 2 Corr. Area 1: 2.3
 Noise: 0.2

Peak table for sample 1 : sample 3

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	48.08	75	48.49	974.3	
3	52.71	121	211.45	2,643.6	
4	113.00	10,380	75.00	10.9	Upper Marker


Region table for sample 1 : sample 3

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	308	18.5	3.33	2.3	1	39.0

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 2 : sample 6

Number of peaks found: 2 Corr. Area 1: 0.8
 Noise: 0.2

Peak table for sample 2 : sample 6

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	48.30	77	16.70	326.8	
3	52.73	121	45.34	565.8	
4	113.00	10,380	75.00	10.9	Upper Marker


Region table for sample 2 : sample 6

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	483	3.4	0.99	0.8	1	22.9

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 3 : sample 11

Number of peaks found: 0 Corr. Area 1: 0.4
 Noise: 0.2

Peak table for sample 3 : sample 11

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	113.00	10,380	75.00	10.9	Upper Marker


Region table for sample 3 : sample 11

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	487	1.9	0.49	0.4	4	37.1

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 4 : sample 13

Number of peaks found: 4 Corr. Area 1: 0.7
 Noise: 0.1

Peak table for sample 4 : sample 13

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	43.82	40	23.90	898.4	
3	47.34	69	20.40	451.1	
4	49.07	85	19.25	345.1	
5	52.71	121	205.03	2,561.2	
6	113.00	10,380	75.00	10.9	Upper Marker


Region table for sample 4 : sample 13

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	296	5.3	0.98	0.7	0	19.6

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 5 : sample 14

Number of peaks found: 1 Corr. Area 1: 0.7
 Noise: 0.2

Peak table for sample 5 : sample 14

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	52.78	122	178.39	2,215.8	
3	113.00	10,380	75.00	10.9	Upper Marker


Region table for sample 5 : sample 14

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	277	5.3	0.91	0.7	0	21.7

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 6 : sample 15

Number of peaks found: 1 Corr. Area 1: 4.8
 Noise: 0.2

Peak table for sample 6 : sample 15

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	52.69	121	57.25	717.0	
3	113.00	10,380	75.00	10.9	Upper Marker


Region table for sample 6 : sample 15

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	360	30.2	6.20	4.8	6	37.3

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 7 : sample 17

Number of peaks found: 2 Corr. Area 1: 0.4
 Noise: 0.1

Peak table for sample 7 : sample 17

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	52.73	121	97.06	1,211.3	
3	113.00	10,380	75.00	10.9	Upper Marker
4	115.33	12,675	0.00	0.0	


Region table for sample 7 : sample 17

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	518	1.6	0.48	0.4	0	28.8

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
 Modified: 2/18/2014 12:01:55 PM

Electropherogram Summary Continued ...



Overall Results for sample 8 : sample 18

Number of peaks found: 5 Corr. Area 1: 2.2
 Noise: 0.1

Peak table for sample 8 : sample 18

Peak	Aligned Migration Time [s]	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	43.00	35	125.00	5,411.3	Lower Marker
2	43.97	41	16.45	604.1	
3	49.11	85	47.55	848.7	
4	52.77	122	65.29	811.7	
5	53.28	127	98.21	1,168.1	
6	113.00	10,380	75.00	10.9	Upper Marker
7	115.41	12,752	0.00	0.0	


Region table for sample 8 : sample 18

From [bp]	To [bp]	Average Size [bp]	Molarity [pmol/l]	Conc. [pg/μl]	Co Corr. lor Area	% of Total	Size distribution in CV [%]
200	1,000	424	11.7	2.60	2.2	1	33.6

Assay Class: High Sensitivity DNA Assay
Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
Modified: 2/18/2014 12:01:55 PM

Gel Image



Assay Class: High Sensitivity DNA Assay
Data Path: C:\...ttings\Bioanalyzer\2014-02-18\2014-02-18\2014-02-18_001.xad

Created: 2/18/2014 11:29:08 AM
Modified: 2/18/2014 12:01:55 PM

Invalid Samples

Sample 9 has not been run, no results available.

Sample 10 has not been run, no results available.

Sample 11 has not been run, no results available.